

Söring

short clinical abstract:

Special field: Surgery

Sub: Thoraxsurgery

Indication: Pleurectomy,
Segmentresection of the
pneuma

Doctor in charge: DR. ENGELMANN, FLT Berlin, Germany

Diagnosis: Diffusiv expanse bleeding in the area of the chest wall (after pleurectomy), as well at the pulmonal surface (after segmentresection of the pneuma)

Use of the Arco™ : The coagulation is contactless, so that a very large area can be treated in a distance from 1-1,5cm like using an air-brush. Especially good for sealing of large parenchymatous haemorrhage, for punctiform seeping haemorrhage and air escape.

Advantage of the Arco™ : The haemostasis is very quick and lasting. The coagulate tissue shows no carbonization, on the contrary it impress with a very thin eschar.

Literatur: Own documents