

Söring

short clinical abstract:

Special field: Gastroenterology

Sub: Gastrointestinal endoscopy

Diagnosis: Stenosing esophagustumor

Doctor in charge: Dr WERNER JOHANNNS, Clinicum Wuppertal, Germany

Conventional Treatment: Injection of epinephrine (adrenaline) and fibrin glue

Use of the Arco™ : The Argon-coagulation (2-3l /min, bzw 40-75W) was used as a palliative measure for reopening the lumen, before implanting an endoprosthesis.

Advantage of the Arco™ : Very good for haemostatic, removal of residual sessile adenoma tissue, for restoring the oesophageal potency.

Used handpieces: Arco®-MC, Argo-flex (2 mm flexible probe)

Literatur: German-Magazine: Zeitschrift für Gastroenterologie, Dezember 1995, Band 33, S. 694 ff, Demeter-Verlag, D- 72336 Balingen